

Hygiene Monitoring System

Instructions for use

Results reference guide

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**NOTE:** Reference colors exhibited in printed prospects and boxes as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

## PRO1 MICRO

Pen system for surface protein and allergen detection

EN

### Product description

Tissue and blood proteins are residual components frequently found on contaminated surgical instruments. Poor cleaning can lead to proteinaceous residues. The residual protein is important because of the continuing risks of transmission of prions (the causative agent of iatrogenic transmissible spongiform encephalopathies such as variant Creutzfeldt-Jakob disease (vCJD)). It is necessary to use protein detection methods to check for the efficient removal of protein from surgical instruments after processing since protein levels are used as an indication of the amount of prion protein contamination. The Hygiene Monitoring System Chemdye® PRO1 MICRO has been designed to detect proteins and to reduce iatrogenic agents on surgical instruments after their cleaning, in real time. The system has a high absorption swab, that allows collection of samples from different surfaces with high efficacy. A visual readout of a color change indicates the presence of detectable levels of protein. Chemdye® PRO1 MICRO has high sensitivity and can detect up to 1.0 µg of protein.

### Characteristics

PRO1 MICRO provides all the necessary equipment required for HTM 01-01/ HTM 01-05 (UK) / ISO 15883-1 Residual Protein Test.

### Indications for use

Use for testing any surface after the cleaning process and for hard-to-reach areas in complex instruments such as serrated edges or box joints.

### Precautions

Do not open until use. Do not reuse. Do not touch the swab with your hands. Wear gloves, surgical mask and other hygiene measures. Do not freeze.

High concentrations of alkaline sanitizers may cause false negative results in some cases. Hydrogen peroxide-based disinfectants may cause a false positive color change reaction. Although Chemdye® PRO1 MICRO has high sensitivity, it is not a microbiological monitoring device, and surfaces may not be considered sterile based on the results obtained.

### Instructions for use

1. Remove the swab from the device (Fig. A). Carefully apply 2 drops of moisturizer on the swab or surface to be tested.

2. On the surface to take the sample, swab thoroughly (Fig. B).

**NOTE:** The surface size to be sampled depends upon the standard that you are following.

A) If you work under the ISO 15883 standard, it is recommended to sample approximately 10 cm<sup>2</sup> of surface. Swab in zigzag in one direction and then do the same perpendicularly (Fig. B). Apply pressure on the swab and rotate it while you collect the sample. For performing comparable cleaning tests, you must standardize surface sampling. In order to achieve this, it is advisable to take the sample always at similar points of surface, with an invariable area of 10 cm<sup>2</sup>, swab 10 times in each direction.

B) If you work under the HTM 01-01 standard, pick one size of the surgical instrument to be tested and swab thoroughly over its surface. Apply pressure on the swab and rotate it while you collect the sample.

3. Place the swab back into the device (Fig. C)

4. Activate it by pressing it down firmly (Fig. D).

5. Slide the swab upwards. (Fig E)

6. Shake vigorously downwards until the solution turns to green. (Fig F)

**CAUTION:** Not shaking the device properly so that the liquid of the inner compartment goes down completely, may bring invalid results.

7. Slide the swab downwards. (Fig. G)

8. Shake vigorously downwards for 15 seconds. (Fig. H)

9. Slide the swab upwards. (Fig. I)

10. Shake vigorously and collect the solution in the readout cone (Fig. J)

11. Make sure that the swab is no longer seen within the readout cone. (Fig. K)

12. Incubate the Chemdye® PRO1 MICRO device at 55-60 °C for:

- 10 minutes in Bionova® MiniPro

- 15 minutes in Bionova® IC10/20, IC10/20FR, IC10/20FRLCD. (Fig. L)

**IMPORTANT:** Incubation must be performed without the swab immersed in the solution.

13. For a qualitative analysis, take the Chemdye® PRO1 MICRO from the incubator and let it cool to room temperature for 5 minutes. The colorimetric result can be interpreted visually by using the table printed on the label of the device. By comparing the color of the test against the table, an estimation of the surface cleanliness can be made:

**OK** - Clean, no further action required.

**Warning!** - Re-clean.

**x/xx** - Re-clean and re-test.

**IMPORTANT:** The estimations must be done within the 5 minutes after taking the device out of the incubator. Results analyzed beyond this period are not valid.

14. For a quantitative analysis, you may use Bionova® IC10/20FR, IC10/20FRLCD or MiniPro Auto-readers. Read the instructions for use of these devices previously.

15. After using the Chemdye® PRO1 MICRO system, reprocess the instrument or the surface used for the monitoring of proteins.

### Storage

Chemdye® PRO1 MICRO system should be stored at temperatures between 2-30 °C in dark places. Shelf life is 24 months when stored refrigerated (4 °C) and 6 months when stored at room temperature (25 °C).

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