
Indicadores Biológicos Biological Indicators

Para la esterilización con VH_2O_2
 For VH_2O_2 sterilization
 Para esterilização por VH_2O_2
 Pour stérilisation avec VH_2O_2
 Per sterilizzazione con VH_2O_2

BT91

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
Parameters determined at time of manufacture according to EN ISO 11138-1: 2006 and IRAM 37102-1: 1999. The showed values are reproducible only under the same conditions under which it were determined.

ISO and USP Compliant
 ATCC is registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
 Director Técnico
 Quality Assurance Director

Industria Argentina - Made in Argentina

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English

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores. It also has an growth indicator media of purple color contained in the glass ampoule.

Product description

Bionova® BT91 biological indicator is specifically designed for the monitoring of Hydrogen Peroxide sterilization processes. If sterilization process was not successful, the indicator media will change from purple to yellow after incubation at 60 ± 2 °C, thus indicating the presence of live *Geobacillus stearothermophilus* spores. If the sterilization process was successful, the indicator media will remain purple after incubation. The final readout should be made after 24 hours of incubation at 60 ± 2 °C.

Warning!

Do not use BT91 biological indicator for monitoring steam sterilization cycle, or ethylene oxide, dry heat, formaldehyde or other sterilization processes.
 Do not re-use the biological indicators.

Storage

Store it in a dark place and temperatures between 15-30 °C, 35-60% relative humidity.
 Do not freeze.
 Do not store these biological indicators near sterilizing agents or other chemical products.

Directions for use

- Identify the Bionova® BT91 indicator by writing the sterilizer number (in case of having more than one sterilizer), load number, and processing date on the indicator label.
- Pack the biological indicator along with materials to be sterilized in an appropriated package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for sterilizing agent.
- Sterilize as usual.
- After the sterilization process has finished, you can:
 - Remove the biological indicator from sterilization package for incubation prior to material aeration.
 - Aerate the package containing the biological indicator with the rest of the load and then remove the biological indicator for incubation.
- Check the chemical indicator on the label of biological indicator. A color change to pink confirms that the biological indicator has been exposed to Hydrogen Peroxide. **IMPORTANT:** this color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.
- Crush the glass ampule contained in the biological indicator and incubate at 60 ± 2 °C.

IMPORTANT: Use a non-sterilized biological indicator as a positive control each time a processed biological indicator is incubated. The positive control ensure that correct incubation conditions were met.

- Incubate the processed biological indicator and the indicator used a positive control for a maximum of 24 hours at 60 ± 2 °C.

A color change from purple to yellow of the growth indicator media means that sterilization process failure has occurred. If after 24 hours there is no color change in the processed indicators, a final negative result is made (the sterilization process was acceptable). The positive control indicator should show a purple to yellow color change for the result to be valid.

NOTE: To extend the incubation time to more than 24 hours, we recommend using a humidified environment to prevent the complete evaporation of the culture medium contained in the indicator.

Record the results and discard immediately as it is shown below.
WARNING! Do not use the sterilizer until the biological indicator test results are negative (process indicator reads purple).

Disposal

Dispose of used biological indicators according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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