KPCD222-2

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KPCD222-2

Steam Test Pack PCD Kit

Process Challenge Device for Steam sterilization processes.

Usage

For routine load release, especially for loads containing implants. For routine sterilizer monitoring. For periodic validation of the sterilizer (after installation, repair, relocation).

Detects inadequate air removal and steam penetration in dynamic-air-removal (pre-vacuum) and gravity displacement steam autoclaves.

Conditions: 132/135 °C ≥ 4 minutes and 121 °C ≥ 30 minutes, steam sterilization processes.

Exclusive conditions for United States:

- Gravity cycles at 121 °C, 30 minutes.
- Gravity cycles at 132 °C, 15/25 minutes.
- Gravity cycles at 135 °C, 10 minutes.
- Vacuum assisted cycles at 132 °C, 4 minutes.
- Vacuum assisted cycles at 135 °C, 3 minutes.

Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/ NS-EN ISO 13485:2016.

ISO 11140-1:2014, ISO 11138-1:2017, ISO 11138-3:2017, ANSI/AAMI ST79.

Classification

Class 1, according to risk.

Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

Characteristics

The kit consist of 25 Porous cards systems plus a Bionova® BT222 Self-Contained Biological Indicator (SCBI) box (25 units). Each porous cards system contains a Bionova® BT222 SCBI and a Type 5 Integrator Indicator printed on PCDBI-2-RC self-adhesive Record Card, where data on sterilization cycle may be written. Two possibilities for documentation record: the whole Record Card can be adhered, or alternatively, only the area containing the indicator may be recorded.

Process Indicator:

Initial color: **blue**.


Final color: **dark grey or black**.

PCDBI-2-RC Integrator Indicator (within the challenge device):

Initial color: **purple**.

Final color: **green**.

100 % Toxic Heavy Metals free.

Integration condition is calibrated against the kill time of a 10^6 *G. stearothermophilus* ATCC 7953 spore population, calculated in BIER (Biological Indicator Evaluator Resistometer). Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

BT222 Self-contained Super Rapid Biological Indicator (within the challenge device): $\geq 10^6$ *G. stearothermophilus* ATCC 7953 spores per vial in paper carrier. Glass ampoule with growth indicator medium. Label with chemical indicator line printed with Steam reactive ink (Color change: pink to brown).

Final results: by fluorescence after incubation for 1 hour at 60 °C (sensitivity $\geq 97\%$).

An optional visual color change confirmation could be made after 48 hours of incubation. A 7-day readout is optional to confirm 1-hour result by color change.

NOTE: If sterilization process was not successful, indicator medium will turn yellow after incubation, thus indicating the presence of living *G. stearothermophilus* spores. If sterilization process was successful, indicator medium will remain purple after incubation.

Environmental conditions during manufacture

15-30 °C (room temperature), RH 30-80 %.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep out of direct light.

Transportation conditions

Storage conditions should be strictly followed. Carry in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

Shelf-life

2 years.

Packing

25 PCD222-2 + 25 BT222 Biological Indicators per box.

Packing information: product description, storage conditions, manufacturer information and data on box's label.

Labelling

On product's box: Product code and description, process for intended use, batch number, batch number of the Self-Contained Biological Indicator, expiration date, normative, barcode and datamatrix.

On the Bionova® BT222's box: product code and description, batch number, load strain, manufacture and expiration date, barcode and datamatrix.

On the Bionova PCD222-2's box: Product's code, process for intended use, batch number, batch number of the Self-Contained Biological Indicator, manufacture and expiration dates, regulation, barcode and datamatrix.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical Industries.

Other important information

Read product's directions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to EO, Dry Heat, Radiation or any sterilization process other than Steam.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement steam sterilization, or at 134 °C for at least 10 minutes in vacuum steam sterilizer.

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