
Results Reference Guide


CDS47V

CDS47A


Unprocessed

Failures

Processed and Satisfactory

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EN Chemical Indicators

For monitoring of Airborne and Surface Disinfection processes with Hydrogen Peroxide Vapor or Aerosol

Product description

Chemical Indicators Chemdye® ChemSurf CDS47A and CDS47V have been designed to react in disinfection processes of airborne and surfaces by Hydrogen Peroxide. CDS47A and CDS47V inks were developed to change when the disinfection process reaches the values indicated for the critical process variables.

Characteristics

Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators consist of a pyramidal synthetic support (10.5 cm on the side and 9.0 cm in height) printed with a specific indicator ink designed to react to aerosol or vapor Hydrogen Peroxide, respectively. Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators present an exclusive design that allows it to dangle in height, according to the need. See forms of use in height and corresponding assemblies.

CDS47V (Vapor)

- Initial color: red
- Final color: yellow

CDS47A (Aerosol)

- Initial color: bordeaux
- Final color: green

Indicators 100% free of toxic heavy metals.

Indications of use

Use Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators to monitor airborne and surface disinfection processes by vapor or aerosol Hydrogen Peroxide.

Conditions for CDS47V (Vapor)*

- H₂O₂ minimum room concentration ≥ 400 ppm.
- Exposure time ≥ 20 minutes.
- Temperature: 10-30 °C, RH: 30-80%.

Conditions for CDS47A (Aerosol)*

- H₂O₂ solution concentration ≥ 6%.

- Injected quantity: 5 ml solution of H₂O₂ / m³ room.
- Exposure time ≥ 2 hours.
- Temperature: 10-30 °C, RH: 30-80%.

(*) By individual enclosure.

Precautions

Do not use Chemdye® ChemSurf CDS47A and CDS47V Chemical Indicators to control disinfection processes with other agents than Hydrogen Peroxide or for sterilization processes by Hydrogen Peroxide in chamber.

Instructions for use

1. Perform the three-dimensional assembly of the indicators, following the assembly instructions specified below.
2. Place the indicators in those areas that you consider a priori more inaccessible to vapor or aerosol hydrogen peroxide. Locations that experience minimal H₂O₂ flows are suggested, such as the corners of the enclosure or behind a bulky item or furniture. To define which the most difficult areas to decontaminate are, it is recommended to use several indicators to validate the same disinfection cycle.
- NOTE: The CSDC documentation chart is attached along with the product. You can indicate there the position of the indicators with respect to the device used for environmental disinfection. Also, Chemdye® ChemSurf CDS47A and CDS47V chemical indicators have a grid on one of its faces, which also allows to locate the position of the indicator and the equipment, individually.
3. Perform the decontamination and airborne process as usual.
4. After the process is finished, analyze the results. The indicators must change to the reference color in all its faces to consider that the adequate disinfection conditions were reached.
- NOTE: If the indicators are placed on a surface, the face in contact with it will not suffer a color change because it won't be exposed to the disinfecting agent.


Storage and shelf life

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store close to sterilizing agents. Indicators have an expiration date of 5 years from the date of manufacture when stored under recommended conditions. Both the expiration date and lot number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

Disposal

Discard chemical indicators after use according to your country's healthcare and safety regulations.

Figure 1

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1. Fold the indicator by the dotted lines, leaving the chemical ink out.
2. Insert the flap "A" in the slot corresponding to the same letter, as indicated in Fig.1. Repeat the same for all the flaps.
3. To hang the indicator, pass the silicone seal through the perforations of faces 2 and 3. Fix the self-adhesive hook to a suitable surface and hang the pyramid from the seal.