KBT400

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KBT400

Spore Strips Biological Indicator Kit

For Liquid Peracetic Acid High Level Disinfection processes.

Usage

Monitoring the sporicidal activity of Liquid Peracetic Acid in High Level Disinfection processes.

Applicable Regulation

ISO 11138-1:2006 and IRAM 37102-1:1999.

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.


Characteristics

25 mm x 70 mm non-absorbent material package. Turquoise code.

5 mm x 25 mm polyethylene fibers strip.

10^6 *Geobacillus stearothermophilus* ATCC 7953 spores per strip.

Silicone fastener: 80 mm long. Pearl end: 5.6 mm diameter. Color: turquoise.

Upon completion of the sterilization cycle, the spore strip should be tested in a suitable culture medium at 60°C. We suggest the use of BIONOVA® MC1020 culture media to final reading in 24 hours. For conventional culture media like TSB, is recommended incubation for 7 days at 60 °C.

D-value: at 50°C, 3.500 ppm of Peracetic Acid.

MC1020-2. Culture medium with pH indicator for growth of *Geobacillus stearothermophilus* ATCC 7953, *Bacillus subtilis* DSM 5230 and *Bacillus atrophaeus* ATCC 9372 bacterial spores.

Growth medium color: red.

Plastic tube: 76.0 mm high, 12.5 mm diameter.

Medium volume: 2 ml.

Medium turns to yellow when spores grow.

Environmental conditions during manufacture

Room temperature 15-30 °C, RH 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

Storage conditions

T = 10-30 °C, RH 30-80 %, keep in a dark place.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

24 months.

Packing

100 spore strips envelopes + 100 culture medium tubes per box.

Weight: 643.8 g.

Packing information: code and description of the product, process for intended use, regulation, manufacturer information and data on box's label.

Information on envelope: code and description of the product, bacterial stain, batch number and expiration date.

Labelling

On product's box: code and description of the product, product batch number, kit components batch number and expiration date.

On foil bag: code and description of the product, regulation, batch number, bacterial population, presentation and manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

It is advisable to incubate at 60 °C in BIONOVA® IC10/20 incubator.

Read product's directions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to radiation, EO or any processes other than Liquid Peracetic Acid.