Bowie & Dick Test Pack

Rev. 9 / 03.2018


Paquete de Prueba Bowie & Dick BD125X/1

ES

Bowie & Dick Test Pack BD125X/1

EN

Unprocessed




Warning Sheet

Main Indicator

Unsatisfactory

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Description

Chemdye® Bowie & Dick Test Pack detects air leaks, inadequate air removal and steam penetration in vacuum-assisted steam sterilizers. The Chemdye® Bowie & Dick Test Pack consists of two chemical indicator sheets (125 x 113 mm) positioned between porous materials and a reticulated foam sheet wrapped all together forming a package, with a steam indicator label.

The central sheet with a regular pattern corresponds to a Bowie & Dick classic indicator, i.e. BD Test Sheet. The chemical indicator changes from yellow to dark brown/black when processed. Any unexpected color change, such as the center of the test sheet being paler or a different color from the edges (i.e., there is a non-uniform color change) indicates that there was an air pocket present during the cycle due to sterilizer malfunction.

The package has also a Warning Sheet that contains a circular chemical indicator. The Warning Sheet is placed at a lower point than the BD Test Sheet inside the package and allows rapid detection of failures in the sterilization cycle.

Bowie & Dick Test Pack should be used daily, before running the first load of the day, after a sterilizer is installed or relocated, after a sterilizer malfunction, after sterilization process failures, and after any major repairs of the sterilizer. A shortened cycle (i.e., a cycle omitting the post-vacuum drying phase) should be run first to properly heat up the sterilizer.

Instructions for use

1. Place the Bowie & Dick Test Pack horizontally with the label facing up in an empty vacuum-assisted sterilizer chamber, in the bottom/front of the sterilizer rack, and over the drain.
2. Run the sterilizer for not more than 3.5 min. - 134 °C / 4 min. - 132 °C.

Interpretation of results

After completing the cycle, retrieve and examine the chemical indicator sheet.

Satisfactory test

It is indicated by a test sheet that shows uniform dark brown/black color development. Such a test indicates rapid steam penetration, adequate air removal and lack of significant air leaks.

Unsatisfactory test

It is caused by an air pocket within the chamber and indicates incomplete air removal from the chamber, poor quality wet steam, superheated steam and/or non-condensable gases in steam supply.

If a non-uniform color change is noted, it should immediately be reported to the supervisor and that sterilizer should not be used until examined for possible malfunction. Recall and quarantine the loads sterilized in that machine until it is certain that the machine has not failed before.

Warning sheet

If there is no evidence of a black or darker color than the reference color in the circle, the equipment must be examined for possible malfunction.

Warning!

Do not use a Bowie & Dick Test Pack in gravity displacement steam sterilization cycles (where there is no pre-vacuum phase to check) or in other kinds of sterilization processes.

Storage

Store Chemdye® Bowie & Dick Test Packs in a dark place, 30-80 % RH at 10-30 °C. It is recommended to protect the indicator Test Sheet after use inside a plastic foil for its storage.

Disposal

Discard along with cardboard waste according to your country's healthcare regulations.

The Bowie & Dick Test Pack is manufactured according to ISO 11140-1:2014 and ISO 11140-5:2007 standards.

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

Innovation in technologies for sterilization and disinfection control

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