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Composición

Cada tubo contains a population of *Bacillus atrophaeus* ATCC 9372 spores soaked on a paper strip (spores carrier). It also contains blue culture medium within a glass ampoule, which is on top of the filter and the carrier, on the base of the tube.

Product description

Bionova® BT110 Rapid Readout Biological Indicators for Ethylene Oxide sterilization processes have been designed for quick and easy monitoring of Ethylene Oxide (EO) sterilization cycles. They consist of a plastic tube, a special filter, a spores carrier and a glass ampoule with culture medium. The tube has a plastic cap with holes and a permeable to EO barrier.

Rapid readout: 4 hours

The Rapid Readout must be carried out in the Bionova® IC10/20FR, Bionova® IC10/20FRLCD or Bionova® MiniBio reader/incubator or in a similar reader/incubator that meets the requirements that are described next. The Reader/incubator must be capable of reading the fluorescence emission of the product, resulting from breaking a specific substrate on the spores carrier. Fluorescence is produced when the reader stimulates the carrier with UV light at 360 nm. Final reading of negative results is readily available at 4 hours of incubation. Fluorescence readout is an indirect determination of the germination and growth of *Bacillus atrophaeus* spores which have survived the sterilization process (positive result). Furthermore, a failure in the sterilization process can also be evident by the color change of the medium. Due to the high sensitivity of the fluorescence results after 4 hours, conventional incubation for color change of BT110 Rapid Readout Indicator is not an advantage.

WARNING: Do not use Bionova® BT110 Biological Indicators to control steam, dry heat, formaldehyde or any other sterilization process than EO. Do not reuse the Biological Indicators.

Storage

Store it in a dark place at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

1. Identify the Bionova® BT110 Biological Indicator by writing the sterilizer number (if you have more than one), load number and processing date on the label.
 2. Pack the Biological Indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which you consider a priori most inaccessible for the sterilizing agent (EO). Generally, a problematic area is the center of the load or near the sterilizer door.
 3. Sterilize as usual.
 4. After the sterilization process has finished, open the sterilizer door and remove the biological indicator from the package.
 5. Check that the chemical Indicator printed on the label of the biological indicator has changed to green. The color change indicates that the biological indicator has been exposed to Ethylene Oxide. **IMPORTANT!** This color change is not an evidence that the process was sufficient to achieve sterility. If chemical Indicator color has not changed, it is necessary to check the sterilization process.
 6. Press the lid to seal the tube. Then, break the ampoule contained within the biological indicator. This may be done in 3 different ways:
 A- With an ampoule crusher.
 B- Manually. To do this, take the tube by placing your index finger and thumb on the space between the cap and the line of chemical Indicator and then press the tube.
 C- With the ampoule crusher placed within the top of the incubator's incubation area.
- Afterwards, shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaked the spores carrier entirely. Finally, place the biological indicator in the incubator.
- Logo, agitate the indicator with similar movements to those used to lower the temperature in a mercury thermometer. Shake the indicator every time a processed Indicator is incubated. The positive control ensures that correct incubation conditions were met; capability of medium to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals; and proper functioning of Bionova® IC10/20FR, Bionova® IC10/20FRLCD or Bionova® MiniBio Indicator. Both the positive control indicator and the processed indicator should belong to the same batch.
7. Incubate the processed biological indicator and the Indicator used as positive control for a maximum of 4 hours at 37 ± 2 °C for rapid readout. Fluorescence detection by the reader (stimulation 340-380 nm/emission 455-465 nm) means failure in the sterilization process. If no fluorescence is detected at 4 hours of incubation (chamber relative humidity ≥ 35%), then the result is negative (the sterilization process has been effective). For results to be valid, the indicator used as positive control must yield positive in the reader. It is good practice to incubate a positive control for a visual color change. Record the results and discard them immediately as it is explained below. **WARNING!** Do not reuse the sterilizer until the biological indicator test result is negative (the processed Indicator medium must remain blue). Once the fluorescence readout of the Biological Indicators has finished, do not restart it, since results might be erroneous. If chamber relative humidity is unknown or < 35%, continue to incubate the Biological indicator for a pH color change result, since low relative humidity may increase the fluorescent readout time beyond 4 hours.

Visual confirmation: 48 hours

Optionally, you can perform a visual color change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will change to a greenish color first, and then to yellow during incubation at 37 ± 2 °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain blue after the incubation process.

A definitive negative result is obtained after the 48-hour incubation.

The positive control must show a color change from blue to yellow for results to be valid.

Visual readout: 7 days

A visual readout may be done when chamber relative humidity is unknown or < 35%, but is optional and not intended to be routinely performed when it is ≥ 35%. This is an initial validation of the 4 hour reading. The fluorescent results are compared to the 7-day visual readings. The sensitivity of the system is determined as the difference between the number of 7-days positive indicators and false negatives indicators (negative for fluorescence reading and visually positive) in relation to the number of 7 days positive Indicators (*). Based on the claimed ≥ 97% readout sensitivity of the 4-hour Biological Indicator, there is no advantage in incubating the indicator beyond 4 hours when the chamber relative humidity is ≥ 35%. **NOTE:** If 7-day readout is performed, a humidified environment will be required to prevent medium from drying.

Disposal

Dispose biological indicators after use, according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, or at 132 °C for 15 minutes in a gravity displacement steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted steam sterilizer.

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