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## Dual Incubator for Biological Indicators and Hygiene Monitoring Systems

### Product description

Bionova® IC10/20 incubator has been designed for a wide range of products. It has a heating block with:  
-34 positions for conventional and Rapid Readout Colorimetric Self-contained Biological Indicators (8 hours).  
-14 positions for Self-contained Ampoules, Culture Medium and Protein Pen Detection System.  
-A special hole for external thermometer.  
-An ampoule crusher for IBs.  
The incubator is pre-set at 37 or 60±2 °C. These two temperatures cannot be modified.

### Operating conditions

Bionova® IC10/20 incubator works at room temperature of 10-30 °C, Relative Humidity of 30-80% and a supply voltage of 12 V Direct Current.

### Power supply specifications:

-Input parameters: 110/220 V, 0,2 A, 50/60 Hz.  
-Output parameters: 12 V, 4 A.

### WARNING

Do not pour any liquid inside.  
Do not immerse into any liquid.  
Only for indoor use.  
Disconnect the power cord before cleaning.  
Do not use abrasive or corrosive cleaners or disinfectants.  
The incubator should be sent to the manufacturer for repair.  
Make sure that the incubator is connected to a properly rated power cord.

### Directions for use

1. Turn on the machine by connecting the end of the AC power supply to the mains and then connect the other end of the source (plug) on the rear of the IC10/20 unit.

**NOTE:** The proper switching-on of the device may be verified by the blinking of the temperature blue light.

2. Select the incubation temperature using the switch button located at the base of the incubator. During heating, the blue light will remain blinking until it reaches the selected temperature.

**IMPORTANT:** Check that the incubation temperature is the correct one for the product to be incubated.

3. Incubate at 37 or 60±2 °C.

Color change of the growth indicating medium (in conventional and Rapid Readout Colorimetric Self-contained Biological Indicators, Self-contained Ampoules or Culture Medium) shows a failure in the sterilization process.

A color change in the Protein Pen Detection System indicates the presence of detectable protein and/or allergens levels.

**NOTE:** For the incubation and detailed analysis of the results, see specifications contained in each product's directions for use.

### Calibration

Bionova® IC10/20 incubator allows to calibrate the incubation temperature with the Bionova® TB-IC1020 thermometer. It is advised that this procedure be carried out only by qualified technicians. To perform the calibration, follow this steps:

1. Turn on the incubator and select the incubation temperature (37 °C or 60 °C). Wait until the selected temperature stabilizes; the blue light will stop blinking.

2. Locate the calibration hole at the bottom of the incubator.

3. Insert a precision flat-head screwdriver of 1,8 mm in the hole.

4. To rise the incubation temperature, slightly rotate the screwdriver in a counterclockwise direction. Wait 2 minutes, after the blue light becomes stable again, and verify the temperature increase.

5. To lower the incubation temperature, rotate the screwdriver in a clockwise direction. Wait 2 minutes, after the blue light becomes stable again, and verify the temperature variation.

**NOTE:** During the calibration process, the light might become intermittent.

**IMPORTANT:** Due to the high sensitive of the calibration hole, it's recommended to verify the temperature increase every time the screwdriver is rotated ¼ of a turn.

### Warranty

Terragene® S.A. guarantees both the quality of the product material components and the quality of its manufacturing process. Should any material or manufacturing faults be detected within the warranty period (1 year from the time of purchase), the only obligation of Terragene® S.A. will be product repairing or substitution.

### Limitation of liability

Terragene® S.A. shall not be liable for any loss or damage that may result from the unsuitable use of the equipment, negligence or user's full responsibility.

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