# BT24 Biological Indicator

Rev. 8 / 11.2020

## Spore Ampoules

BT24/4 - BT24/5 - BT24/6

Fuchsia/Red = Sterile

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### Quality certification Certificado de calidad Bionova® BT24/

Steam sterilization at low temperature /  
Esterilización por Vapor a baja temperatura

*Bacillus subtilis* DSM 5230 ATCC 35021

LOT

Population / Población \_\_\_\_\_ CFU / UFC

D - value / Valor D (121 °C) \_\_\_\_\_ min.

Survival time / Tiempo de sobrevivida \_\_\_\_\_ min.  
Survival time =  $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte \_\_\_\_\_ min.  
Kill time =  $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

D - value / Valor D (115 °C) \_\_\_\_\_ min.

Survival time / Tiempo de sobrevivida \_\_\_\_\_ min.  
Survival time =  $(\log_{10} \text{ labeled population} - 2) \times \text{labeled D-value}$

Kill time / Tiempo de muerte \_\_\_\_\_ min.  
Kill time =  $(\log_{10} \text{ labeled population} + 4) \times \text{labeled D-value}$

Z-value / Valor Z \_\_\_\_\_ °C

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-3:2017 and IRAM 37102:1999 (Parts 1 and 3). The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant.  
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto  
Director Técnico  
Technical Director

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## EN Biological Indicators

For Steam sterilization at low temperature

**Composition**  
Bionova® BT24 Biological Indicator consists of a glass ampoule containing *Bacillus subtilis* DSM 5230 ATCC 35021 spore population dispersed in a culture medium especially designed to change its color when microbial growth occurs.

**Product description**  
Bionova® BT24 Biological Indicator is specifically designed for monitoring of steam sterilization processes at 110 -121 °C.  
If the sterilization process was not successful, the indicator medium will turn to yellow after incubation between 35-39 °C, thus indicating the presence of living *Bacillus subtilis* spores.  
If the sterilization process was effective, the indicator medium will remain the original color after incubation. The final readout should be performed after 48 hours of incubation between 35-39 °C.

**Precautions**  
Do not use Bionova® BT24 Biological Indicator for monitoring Ethylene Oxide, Dry Heat, Chemical Steam, Radiation or other sterilization processes different from its intended use.  
Do not reuse biological indicators. Must be discarded after use.  
The ampoules can be left out for a maximum of two hours at room temperature below 25 °C.

Protect from shocks, falls and sudden temperature changes.  
Check the color of the indicator before starting a sterilization process. It must be red. Discard it otherwise. The negative control vial can be used to compare the color.  
Do not use an indicator after the expiration date indicated in the package.  
The indicator must be processed by a trained healthcare professional.  
Do not use the sterilizer until the biological indicator test results are negative (process indicator remains the original color).

**Storage**  
Store it in a dark place and refrigerated at temperatures between 2-8 °C, 30-80 % relative humidity.  
Do not freeze.  
Do not store these biological indicators near sterilizing agents or other chemical products.

**Instructions for use**  
1. Take the Bionova® BT24 Ampoules out of the refrigerator and leave them at room temperature during not longer than ten minutes.  
2. Place the Bionova® BT24 Ampoules inside a container with the liquid being sterilized or in a container with liquid similar to that containing the liquid to be sterilized. This container should be of the same size and have the same volume as the one containing the liquid being sterilized.  
3. Sterilize as usual.  
4. Once the sterilization process has finished, open the sterilizer door, wait five minutes and then remove the biological indicator from the container. **WARNING:** Wear safety glasses and gloves when removing the Bionova® BT24 Biological Indicator from the sterilized container.  
5. Incubate the processed biological indicator, a negative control and a non-sterilized biological indicator as a positive control, during 48 hours between 35-39 °C.  
The positive control ensures that the incubation conditions were met.  
The negative control, not being exposed to the sterilization process, can be used as a reference color to compare positive and negative results. Conversely, when used with the ampoules during the process, it ensures that any color change produced immediately after sterilization is the result of the process and not due to germination and/or spore growth.  
Readout should be performed at convenient intervals of 10 hours.  
It is recommended the use of Bionova® IC10/20 Dual Incubator or similar devices that assure an effective temperature control. **WARNING:** Traditional culture stoves have thermal profiles that differ from the above recommended incubators, thus their use might affect or delay the Bionova® BT21 Biological Indicators result.  
A color change to yellow of the processed biological indicator medium means a sterilization process failure has occurred. If after 48 hours there is no color change in the processed indicators, a final negative result is obtained, meaning the sterilization process was successful. The positive control indicator should turn to yellow for the result to be valid and the negative control should remain the same color.  
Record the positive results and discard them immediately as it is shown below.

**Disposal**  
Discard Bionova® BT24 Biological Indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be sterilized in gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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