Two level Chemical Integrator

Quality Certification

EO

LOT

MW

Exp

Integration conditions

Gas concentration _____

Relative Humidity _____

Stated Value (54 °C) - SAFE condition _____ min.

Stated Value (37 °C) - SAFE condition _____ min.


* Integration condition may be extrapolated to the kill time of *Bionova*® Biological Indicator for EO.

According to ISO 11140-1:2014 standard

Lic. Adrián J. Rovetto
Technical Director

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SAFE condition

Explanation of Symbols

EO Product designed for use with Ethylene Oxide sterilization cycles.

CE CE mark.

EC REP Authorized representative in the European Community.

LOT

Batch number.

MW

Manufacture Date.

Exp

Expiration Date.

Man

Manufacturer.

Intended Use Table

| Ethylene Oxide | |
|----------------|--|
| Model | Cycles |
| IT12 | 55 °C for 240 minutes, 600 mg/l, RH 60 % |
| | Minimum Stated Values (SV) for Chemical Integrators |
| | 37 °C 54 °C 55 °C |
| | 75 minutes 30 minutes 28 minutes |

Terragene® and Integron® are registered trademarks of Terragene S.A.

EN Two level Chemical Integrator For Ethylene Oxide sterilization processes

Indications for Use United States

The integrator Terragene® Integron® IT12 is designed to chemically react over time with the critical parameters of ethylene oxide sterilization cycle within a specified tolerance. The integrating indicator strip is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles: 55 °C for 240 minutes, 600 mg/l EO, RH 60%. SV37 °C/75 min., SV54 °C/30 min., SV55 °C/28 min. See Intended Use Table for reference.

Outside the United States

Integron® IT12 two level chemical integrator indicator is designed to monitoring Ethylene Oxide sterilization processes.

Device Description

Integron® IT12 Chemical Integrators (Type 5 according to ISO 11140-1:2014 standard) are single-use two level chemical indicators that consist of paper strips laminated on the back side and printed with indicator ink. They have been designed for monitoring Ethylene Oxide (EO) sterilization processes. The integration level contains a purple/brown ink which was developed to turn to green when a theoretical spore population reaches its kill time, thus indicating that the SAFE condition has been met. Level 1 - Red: Exposure level. Indicator ink turns yellow to indicate exposure to the process, but does not ensure that the variables required to achieve the integration condition have been achieved (UNSAFE condition). Level 2 - Purple / brown: Integration level. Indicator ink turns green to indicate that the integration condition indicated in the quality certificate has been achieved. This condition consists of the kill time of a theoretical 10⁶ *Bacillus atrophaeus* spore population, calculated in a BIER (Biological Indicator Evaluator Resistometer), and is called SAFE condition.

Indicators 100 % free of toxic heavy metals

Precautions

Do not use Integron® IT12 Chemical Integrator to monitor Dry Heat, Steam or sterilization methods other than EO. Do not reuse. Do not reuse the sterilizer until the integrator color reaches SAFE condition.

Instructions for Use

1. Place a Integron® IT12 Chemical Integrator along with the material to be sterilized in a suitable package according to the recommended sterilization techniques. Place the indicators in those areas of most difficult access to the sterilizing agent.
2. Perform the sterilization cycle.
3. After the sterilization process has finished and evacuation/ventilation process is over, you may:
A. Remove the integrator indicator taking the necessary precaution to avoid contact with the remaining sterilizing agent.
B. Ventilate the load together with the integrator indicator and then perform the reading.
4. Analyze the results. If the SAFE integration condition has been reached, the sterilization process has been effective. For reference color, please refer to Result Reference Guide.
NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate chemical indicator.

Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

Shelf Life

Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.
Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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