# BT70

## Spore Strip Biological Indicator

For Gamma radiation.

### Usage

Monitoring Gamma Radiation Sterilization Processes.

### Applicable Regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017; IRAM 37102-1:1999.

### Classification

Class 1, according to risk.

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Characteristics

25 mm x 70 mm non-absorbent material package. Fuchsia code.

5 mm x 25 mm filter paper strip.

$10^5$ ,  $10^6$  or  $10^7$  *Bacillus pumilus* ATCC 27142 spores per strip.

Upon completion of the sterilization cycle, the spore strip should be tested in a suitable culture medium at 37 °C. We suggest the use of BIONOVA® MC1020-2 culture medium to final reading in 48 hours. For conventional culture media like TSB, incubation for 7 days at 37 °C is recommended.

D-Value is provided. Calculated using Limited Holcomb-Spearman-Karber method.

### Environmental conditions during manufacture

T= 15-30 °C, RH= 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

T= 10-30 °C, RH= 30-80 %.

### Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

### Shelf life

2 years.

### Packing

100 units per aluminium foil bag.

### Labelling

On product's envelope: product code and description, process for intended use, batch number, expiration date, bacterial strain and datamatrix code.

On product's packing: product code and description, process for intended use, strain and population of bacteria, batch number, manufacture and expiration date, presentation, regulation, storage conditions, manufacturer information, barcode and datamatrix code.

### Possible target markets

Healthcare and Industry.

### Other important information

It is advisable to incubate at 37 °C in BIONOVA® IC10/20 incubator.  
Read product's instructions for use thoroughly before use.

### Precautions

Do not store the product near sterilizing agents.

Do not expose this product to Steam, H<sub>2</sub>O<sub>2</sub> or any sterilization process other than Gamma Radiation.

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Valor D es provisto. Calculado usando el método Limited Holcomb-Spearman-Karber.