BT96

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BT96

Self-Contained Biological Indicator.

Super Rapid Readout Fluorescence System.

Usage

Monitoring plasma or vaporized Hydrogen Peroxide sterilization processes (H₂O₂).

Applicable Regulation

ISO 11138-1:2006 and IRAM 37102-1:1999.

Classification

Class 1, according to risk.

Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012. ANMAT (Argentinean National Administration of Drugs, Food and Medical Device) PM 1614-1.

Characteristics

Polypropylene tube: 50.5 mm high x 8.5 mm external diameter.

Wall thickness: 0.5 mm

Special filter: 17mm of external diameter polypropylene microfibers.

Polypropylene cap: 10.5 mm external diameter, 16.5 mm high. Wall thickness: 1.0 mm

Glass ampoule: 3.8 – 40.0 mm high. mm. External diameter: 6.5 mm.

Wall thickness: 0.2 – 0.3 mm.

Culture medium 0.65 – 0.75 ml, purple color.

Cap filter: polyethylene fibers, 160 g/m², 16.0 mm diameter.

Polyethylene fibers: 5.0 x 40.0 mm, 160 g/m².

≥ 10⁶ *Geobacillus stearothermophilus* ATCC 7953 spores per vial.

Final fluorescence reading is performed after 30-minute incubation at 60 °C (sensitivity: 97 %).

An optional visual pH color change confirmation could be performed after 48 hours of incubation. If sterilization process has not been successful, culture medium will change to yellow during incubation at 60 °C, thus showing the presence of living spores. If sterilization process is successful culture medium will remain purple after the incubation process.

7-day readout is optional and not intended to be routinely performed; it is an initial validation of the 30-minute reading. Fluorescence results may be compared to the 7-day visual reading.

NOTE: If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

D-Value at 2 mg/l H₂O₂, 50 °C.

Environmental conditions during manufacture

Room temperature 15-30°C, relative humidity 30-80%. Sterility conditions are only necessary during the inoculation and manufacturing process, which is performed in laminar flow.

Storage conditions

Best stored in original box under the next conditions: temperature 10-30 °C, relative humidity 30-80 %.

Transportation conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

Shelf-life

24 months.

Packing

10 units per aluminium foil bag. Weight: 35,0 g

50 units per box. Weight: 145,8 g.

Packing information: product code and description, process for intended use, regulation, manufacturer information and data on pack's label.

NOTE: manufacture date is calculated by subtracting 24 months to the expiration date.

Labelling

On product: 17.0 mm x 33.0 mm polypropylene label. 1.5 mm chemical indicator line, printed with Hydrogen Peroxide reactive ink (color change to green). Graph showing final fluorescence reading time, product code, batch number, expiration date, process for intended use and bacterial strain printed in black.

On product's pack: product code and description, batch number, bacterial strain and population, presentation, manufacture and expiration date.

Possible target markets

Healthcare and Industry.

Other important information (optional)

It is advisable to incubate at 60 °C in Bionova® IC10/20FR, IC10/20FR LCD o Mini Bio Auto-Reader incubators.

Read product's directions for use thoroughly before use.

Description of identified and/or predictable risks which could represent a risk for the product

Do not store the product near sterilizing agents.

Do not expose this product to Ethylene Oxide, Dry heat, Radiation or any sterilization process other than Hydrogen Peroxide.

The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes, at 132 °C for at least 15 minutes in gravity displacement Steam sterilization, or at 134 °C for at least 10 minutes in vacuum Steam sterilizer

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